**WFUMB Position Statement:**  
**Equipment cleaning and safe performance of ultrasound examinations**  
**in the**  
**context of COVID-19**

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Abramowicz JS\*<sub>1</sub> and Basseal JM\*<sub>2</sub> (\*co-authors)

<sub>1</sub> WFUMB and Department of Obstetrics and Gynecology, University of Chicago, Chicago, USA

<sub>2</sub> Discipline of Infectious Diseases and Immunology, Faculty of Medicine and Health, The University of Sydney, NSW Australia and Australasian Society for Ultrasound in Medicine, Sydney, NSW Australia

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Canon/Toshiba	<a href="https://global.medical.canon/products/ultrasound/more_information/guideforcleaning">https://global.medical.canon/products/ultrasound/more_information/guideforcleaning</a>
GE	<a href="https://www.gehealthcare.com/products/ultrasound/ultrasound-transducers">https://www.gehealthcare.com/products/ultrasound/ultrasound-transducers</a>
Mindray	<a href="https://www.mindraynorthamerica.com/wp-content/uploads/2019/03/Mindray_M-transducers_disinfection-guide_40369A.pdf">https://www.mindraynorthamerica.com/wp-content/uploads/2019/03/Mindray_M-transducers_disinfection-guide_40369A.pdf</a>
Philips	<a href="https://www.usa.philips.com/c-dam/b2bhc/master/whitepapers/ultrasound-care-and-cleaning/disinfectant-tables-manuals/dt-us.pdf">https://www.usa.philips.com/c-dam/b2bhc/master/whitepapers/ultrasound-care-and-cleaning/disinfectant-tables-manuals/dt-us.pdf</a>
Samsung	<a href="https://samsunghealthcare.com/en/products/uss/RS80A_with_Prestige/Radiology/transducers">https://samsunghealthcare.com/en/products/uss/RS80A with Prestige/Radiology/transducers</a>
Siemens	<a href="https://www.siemens-healthineers.com/en-us/ultrasound/ultrasound-transducer-catalog#Care">https://www.siemens-healthineers.com/en-us/ultrasound/ultrasound-transducer-catalog#Care</a>
Zonare/Mindray	<a href="https://www.mindraynorthamerica.com/wp-content/uploads/2019/03/Mindray_Zonare-transducers_disinfection-guide_40368A.pdf">https://www.mindraynorthamerica.com/wp-content/uploads/2019/03/Mindray_Zonare-transducers_disinfection-guide_40368A.pdf</a>

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- Aakjær Andersen C, Holden S, Vela J, Skovdal Rathleff M, Bach Jensen M. Point-of-Care Ultrasound in General Practice: A Systematic Review. *Ann Fam Med* 2019;12:61-69.
- Abramowicz JS, Evans DH, Fowlkes JB, Marsal K, terHaar G, Committee WS. Guidelines for Cleaning Transvaginal Ultrasound Transducers Between Patients. *Ultrasound Med Biol* 2017;43:1076-79.
- ACEP. American College of Emergency Physicians Ultrasound Guidelines: Emergency Ultrasound Guidelines. *Ann Emerg Med* 2009;53:550-70.
- ACEP. American College of Emergency Physicians (ACEP) Guideline for Ultrasound Transducer Cleaning and Disinfection. *Ann Emerg Med* 2018;72:45-47.
- AIUM. American Institute of Ultrasound in Medicine (AIUM) practice guideline for the performance of obstetric ultrasound examinations. *J Ultrasound Med* 2018;37:1587-96.
- ASUM. Australasian Society for Ultrasound in Medicine Standards of Practice. <https://www.asum.com.au/standards-of-practice/> 2014-2020.
- Bai Y, Yao L, Wei T, Tian F, Jin DY, Chen L, Wang M. Presumed Asymptomatic Carrier Transmission of COVID-19. *JAMA* 2020;10.1001/jama.2020.2565:e202565.
- Basseal J, Westerway S, Juraja M, van de Mortel T, McAuley T, Rippey J, 30-40. eaAJUM. Guidelines for Reprocessing Ultrasound Transducers. . *Australas J of Ultrasound in Med* 2017;20:30-40.
- Boelig R, Saccone G, Bellussi F, Berghella V. MFM Guidance for COVID-19. *Am J Obstet Gynecol MFM* 2020;<https://doi.org/10.1016/j.ajogmf.2020.100106>.
- CDC. 2015 Immediate Need for Healthcare Facilities to Review Procedures for Cleaning, Disinfecting, and Sterilizing Reusable Medical Devices. CDC Health Alert Network. <http://emergency.cdc.gov/han/han00382.asp>.
- Chan KW, Wong VT, Tang SCW. COVID-19: An Update on the Epidemiological, Clinical, Preventive and Therapeutic Evidence and Guidelines of Integrative Chinese-Western Medicine for the Management of 2019 Novel Coronavirus Disease. *Am J Chin Med* 2020:1-26.

- Dietrich C, Goudie A, Chiorean L, Cui X, Gilja O, Dong Y, Abramowicz J, Vinayak S, Westerway S, Nolsøe C, Chou Y, Blaiwas M. Point of Care Ultrasound: A WFUMB Position Paper. *Ultrasound Med Biol*. 2017;43:49-58.
- Gorbalenya A, Baker S, Baric R, de Groot R, Drosten C, Gulyaeva A, Haagmans B, Lauber C, Leontovich A, Neuman B, Penzar D, Perlman S, Poon L, Samborskiy D, Sidorov I, Sola I, Ziebuhr J. Coronaviridae Study Group of the International Committee on Taxonomy of Viruses: The species Severe acute respiratory syndrome-related coronavirus: classifying 2019-nCoV and naming it SARS-CoV-2. *Nature microbiology* 2020;10.1038/s41564-020-0695-z. Advance online publication. <https://doi.org/10.1038/s41564-020-0695-z>.
- Holshue ML, DeBolt C, Lindquist S, Lofy KH, Wiesman J, Bruce H, Spitters C, Ericson K, Wilkerson S, Tural A, Diaz G, Cohn A, Fox L, Patel A, Gerber SI, Kim L, Tong S, Lu X, Lindstrom S, Pallansch MA, Weldon WC, Biggs HM, Uyeki TM, Pillai SK. First Case of 2019 Novel Coronavirus in the United States. *N Engl J Med* 2020;382:929-36.
- Hoyer R, Adhikari S, Amini R. Ultrasound transducer disinfection in emergency medicine practice. *Antimicrob Resist Infect Control* 2016;5:12.
- Kampf G, Todt D, Pfaender S, Steinmann E. Persistence of coronaviruses on inanimate surfaces and their inactivation with biocidal agents. *J Hosp Infect* 2020;104:246-51.
- Kanagala P, Bradley C, Hoffman P, Steeds RP, British Society of E. Guidelines for transoesophageal echocardiographic probe cleaning and disinfection from the British Society of Echocardiography. *European journal of echocardiography : the journal of the Working Group on Echocardiography of the European Society of Cardiology* 2011;12:i17-23.
- Lichtenstein DA. Lung ultrasound in the critically ill. *Ann Intensive Care* 2014;4:1.
- Liu R, Nomura J, Tayal V. POLICY STATEMENT Guideline for Ultrasound Transducer Cleaning and Disinfection. *Annals of Emergency Medicine* 2018;72:1016.
- Muller T, Martiny H, Merz E, Doffert J, Wustner M, Lessel W, Heynemann H, Enzmann T, Dudwiesus H, Nuernberg D, Tesch C, Weber MA, Krishnabhakdi S, Heil J, Wree A, Jenssen C. DEGUM Recommendations on Infection Prevention in Ultrasound and Endoscopic Ultrasound. *Ultraschall Med* 2018;39:284-303.
- Nyhsen CM, Humphreys H, Koerner RJ, Grenier N, Brady A, Sidhu P, Nicolau C, Mostbeck G, D'Onofrio M, Gangi A, Claudon M. Infection prevention and control in ultrasound - best practice recommendations from the European Society of Radiology Ultrasound Working Group. *Insights Imaging* 2017;8:523-35.
- Peng Q, Wang X, Zhang L. Findings of lung ultrasonography of novel corona virus pneumonia during the 2019-2020 epidemic. *Intensive Care Med* 2020;doi:10.1007/s00134-020-05996-6.
- Poggiali E, Dacrema A, Bastoni D, Tinelli V, Demichele E, Mateo Ramos P, Marciano T, Silva M, Vercelli A, Magnacavallo A. Can Lung US Help Critical Care Clinicians in the Early Diagnosis of Novel Coronavirus (COVID-19) Pneumonia? *Radiology* 2020;10.1148/radiol.2020200847:200847.
- ScOR/BMUS. 2019 Society and College of Radiographers and British Medical Ultrasound Society Guidelines For Professional Ultrasound Practice. 4th ed. edn. London, UK, 1-146.
- Skowronek P, Wojciechowski A, Leszczynski P, Olszewski P, Sibinski M, Polgaj M, Synder M. Can diagnostic ultrasound scanners be a potential vector of opportunistic bacterial infection? *Med Ultrason* 2016;18:326-31.
- Soldati G, Smargiassi A, Inchingolo R, Buonsenso D, Perrone T, Briganti DF, Perlini S, Torri E, Mariani A, Mossolani EE, Tursi F, Mento F, Demi L. Is there a role for lung ultrasound during the COVID-19 pandemic? *J Ultrasound Med* 2020;<https://doi.org/10.1002/jum.15284>.
- Westerway SC, Basseal JM. The ultrasound unit and infection control - Are we on the right track? *Ultrasound* 2017;25:53-57.

Westerway SC, Basseal JM, Abramowicz JS. Medical Ultrasound Disinfection and Hygiene Practices: WFUMB Global Survey Results. *Ultrasound Med Biol* 2019;45:344-52.  
Yuen KS, Ye ZW, Fung SY, Chan CP, Jin DY. SARS-CoV-2 and COVID-19: The most important research questions. *Cell Biosci* 2020;10:40.

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